

## JB-1A Magnetic Stirrer



### FEATURES

Simple and interesting magnetic stirrer, using an electrical motor spins the magnet module, stirring speed from 0 to 1250rpm.

### SPECIFICATION

| Model              | JB-1A             |
|--------------------|-------------------|
| Stirring Volume    | 0~2000mL          |
| Stirring Speed     | 0~1250rpm         |
| Top Plate Size     | Dia.145mm         |
| Material           | PC                |
| Stir Bar           | 30(L)×7(Dia.)mm   |
| Power Requirements | AC 220V/50Hz      |
| Dimensions         | 185(Dia.)×75(H)mm |
| Weight             | 0.6kg             |

## MS series Magnetic Stirrer



### FEATURES

- High performance hot plate magnetic stirrer comes with a temperature probe.
- Large LCD display clearly shows the timer, temperature and running status.
- Push-button speed control effectively avoids the hot-liquids hurt the operator.
- Adjustable heating temperature, stirring times and speeds.
- Automatic constant temperature through an appendant temperature probe.
- 1 point push-button calibration improves the accuracy of temperature control.

### STIRRER INCLUDES

- Stir bar
- Temperature probe

### SPECIFICATION

| Model              | MS200                  | MS300   | MS400   |
|--------------------|------------------------|---------|---------|
| Stirring Volume    | 0~2000mL               |         |         |
| Stirring Speed     | 0~1250rpm              |         |         |
| Temperature Range  | —                      | 0~300°C | 0~400°C |
| Timer Range        | 0 to 999 minutes       |         |         |
| Top Plate Size     | 135(L)×135(W)mm        |         |         |
| Material           | Stainless Steel        |         |         |
| Stir Bar           | 30(L)×7(Dia.)mm        |         |         |
| Power Requirements | AC 220V/50Hz           |         |         |
| Dimensions         | 230(L)×180(W)×120(H)mm |         |         |
| Weight             | 2.2kg                  |         |         |