

JB-1A Magnetic Stirrer



Features:

JB-1A cartoon magnetic stirrer is a fascinating small instrument which used for stirring the non-viscous liquids. The stirrer uses an electrical motor spinning the permanent magnet, its speeds from 0 to 1250 rpm and can quickly stir liquid in the short time.

Instrument Includes:

Stir bar

Specifications:

Model	JB-1A
Stirring Volume	0~2000mL
Stirring Speed	0~1250rpm
Top Plate	Dia.145mm
Material	PC
Stir Bar	30(L)×7(Dia.)mm
Power Requirements	AC 220V/50Hz
Dimensions	185(Dia.)×75(H)mm
Weight	0.6Kg

MS200/300/400 Magnetic Stirrer



Features:

- High performance hot plate magnetic stirrer comes with a temperature probe.
- Large LCD display shows the timer, temperature and running status indicators.
- Push-button speed control effective avoids the hot-liquids hurt operator.
- Adjustable heating temperature, stirring times and stirring speeds.
- Automatic constant temperature through an appendant temperature probe.
- 1 point push-button calibration improves the accuracy of temperature control.

Instrument Includes:

- Stir bar
- Temperature probe

Specifications:

Model	MS200	MS300	MS400
Stirring Volume	0~2000mL		
Stirring Speed	0~1250rpm		
Temperature Range	—	300°C	400°C
Timer	0~999 minutes		
Top Plate	135(L)×135(W)mm		
Material	Stainless Steel		
Stir Bar	30(L)×7(Dia.)mm		
Power Requirements	AC 220V/50Hz		
Dimensions	230(L)×180(W)×120(H)mm		
Weight	2.2Kg		