

## JB-1A Magnetic Stirrer



### FEATURES

Simple and interesting magnetic stirrer, using an electrical motor spins the magnet module, stirring speed from 0 to 1250rpm.

### SPECIFICATION

Model	JB-1A
Stirring Volume	0~2000mL
Stirring Speed	0~1250rpm
Top Plate Size	Dia.145mm
Material	PC
Stir Bar	30(L)×7(Dia.)mm
Power Requirements	AC 220V/50Hz
Dimensions	185(Dia.)×75(H)mm
Weight	0.6kg

## MS series Magnetic Stirrer



### FEATURES

- High performance hot plate magnetic stirrer comes with a temperature probe.
- Large LCD display clearly shows the timer, temperature and running status.
- Push-button speed control effectively avoids the hot-liquids hurt the operator.
- Adjustable heating temperature, stirring times and speeds.
- Automatic constant temperature through an appendant temperature probe.
- 1 point push-button calibration improves the accuracy of temperature control.

### STIRRER INCLUDES

- Stir bar
- Temperature probe

### SPECIFICATION

Model	MS200	MS300	MS400
Stirring Volume	0~2000mL		
Stirring Speed	0~1250rpm		
Temperature Range	—	0~300°C	0~400°C
Timer Range	0 to 999 minutes		
Top Plate Size	135(L)×135(W)mm		
Material	Stainless Steel		
Stir Bar	30(L)×7(Dia.)mm		
Power Requirements	AC 220V/50Hz		
Dimensions	230(L)×180(W)×120(H)mm		
Weight	2.2kg		