

JB-1A Mini Magnetic Stirrer



Features

Simple and interesting mini magnetic stirrer, using an electrical motor spins the magnet modules, stirring speed from 0 to 1250rpm.

Specifications

Model	JB-1A
Stirring Volume	0~2000mL
Stirring Speed	0~1250rpm
Top Plate Size	Dia.145mm
Material	PC
Stir Bar Size	30(L)×7 (Dia.)mm
Power Requirements	AC 220V/50Hz
Dimensions	185(Dia.)×75(H)mm
Weight	0.6kg

MS Series Hotplate Magnetic Stirrer



Features

- High-performance hotplate magnetic stirrer comes with a temperature probe
- Large LCD display clearly shows the timer, temperature and running status
- Automatic constant temperature through a connected sensor
- 1 point offset calibration ensures the accurate temperature control
- Push-button speed control effectively avoids that hot-liquids hurt the operator
- Adjustable heating temperature, stirring times and speeds

Ordering Information

- MS200: Stirrer and stir bar
- MS300/400: Stirrer, stir bar and temperature probe

Specifications

Model	MS200	MS300	MS300
Stirring Volume	0~2000mL		
Stirring Speed	0~1250rpm		
Heater	—	380W	450W
Hotplate Temperature	—	Max. 300°C	Max. 400°C
Timer Range	Up to 999 minutes		
Top Plate Size	135(L)×135(W)mm		
Material	Stainless steel		
Display	Custom LCD (95×35 mm)		
Power Requirements	AC 220V/50Hz		
Ambient Temperature	< 50°C		
Dimensions	230(L)×180(W)×120(H)mm		
Weight	2.2kg		